

Assignment 2 (Graphics and Animation)

Tasks under this assignment –

- First of all, I downloaded the CT dataset [1] provided by professor Kumar in Lecture 08 (Volume Rendering) that comprised of DICOM images of a heart.
- Next, I began with the code provided in the same lecture to read the DICOM image into vtk using the vtkDICOMReader to read the images and display it using spyder IDE.
- Second task was to change the color transfer function as well as the opacity transfer function under the code as predefined.
- Furthermore, I divided the display into three viewports by providing the x_min, y_min, x_max and y_max values for the display. On the first viewport, I displayed the volume rendering data of the CT dataset. Next, on the second viewport, I displayed the mhd file of the heart provided under Lecture 05 [2] given by professor Kumar on eclass. Finally, on the last viewport, I displayed a combination of both volume rendered data as well as iso surface rendered data.
- Finally, I synced all three viewports using the SetActiveCamera() function of vtk.

1. Information about the **CT dataset**:

The CT dataset comprises of **188 DICOM images** of the heart.

The dimensions of the DICOM images under CT dataset are – **(512 * 512 * 188)**.

The size of the dataset is **49283072 bytes**.

The **voxel resolution** is denoted by the following array –

(x_min, x_max)	(0, 249.512)
(y_min, y_max)	(0, 249.512)
(z_min, z_max)	(0, 130.899)

2. Source Code

```
1 import vtk
2
3 # path to the stored images
4 path = r'CT'
5
6 # Read data
7 reader = vtk.vtkDICOMImageReader()
8 reader.SetDirectoryName(path)
9 reader.Update()
10
11 # Create colour transfer function
12 colorFunc = vtk.vtkColorTransferFunction()
13 colorFunc.AddRGBPoint(-1024, 0.0, 0.0, 0.0)
14 colorFunc.AddRGBPoint(-77, 0.5, 0.2, 0.1)
15 colorFunc.AddRGBPoint(100, 0.9, 0.6, 0.3)
16 colorFunc.AddRGBPoint(180, 1, 0.8, 0.9)
17 colorFunc.AddRGBPoint(260, 0.6, 0.1, 0)
18 colorFunc.AddRGBPoint(3071, 0.7, 0.8, 1)
19
20 # Create opacity transfer function
21 alphaChannelFunc = vtk.vtkPiecewiseFunction()
22 alphaChannelFunc.AddPoint(-1024, 0.0)
23 alphaChannelFunc.AddPoint(-77, 0.0)
24 alphaChannelFunc.AddPoint(179, 0.1)
25 alphaChannelFunc.AddPoint(260, 0.4)
26 alphaChannelFunc.AddPoint(3071, 0.7)
27
28 # View 1 for volume
```

Figure 1 Reading the data, CT function, Opacity function

```
28 # View 1 for volume
29 volume1 = vtk.vtkVolume()
30 ren1 = vtk.vtkRenderer()
31 txt1 = vtk.vtkTextActor()
32 txt1.SetInput("Volume rendering!")
33 txtprop=txt1.GetTextProperty()
34 txtprop.SetFontFamilyToArial()
35 txtprop.SetFontSize(18)
36 txtprop.SetColor(1,1,1)
37
38 # Define volume mapper
39 volumeMapper1 = vtk.vtkSmartVolumeMapper()
40 volumeMapper1.SetInputConnection(reader.GetOutputPort())
41
42 # Define volume properties
43 volumeProperty1 = vtk.vtkVolumeProperty()
44 volumeProperty1.SetScalarOpacity(alphaChannelFunc)
45 volumeProperty1.SetColor(colorFunc)
46 volumeProperty1.ShadeOn()
47
48 # Set the mapper and volume properties
49 volume1.SetMapper(volumeMapper1)
50 volume1.SetProperty(volumeProperty1)
51
52 # create volume renderer
53 ren1.AddVolume(volume1)
54 ren1.AddActor(txt1)
55
```

Figure 2 Volume rendering

```

58# View 2 for ISO mapper
59ren2 = vtk.vtkRenderer()
60txt2 = vtk.vtkTextActor()
61txt2.SetInput("ISO rendering!")
62txtprop=txt2.GetTextProperty()
63txtprop.SetFontFamilyToArial()
64txtprop.SetFontSize(18)
65txtprop.SetColor(1,1,1)
66
67# Apply Marching Cubes algorithm
68march_cubes = vtk.vtkMarchingCubes()
69march_cubes.SetInputConnection(reader.GetOutputPort())
70march_cubes.ComputeNormalsOn()
71march_cubes.ComputeGradientsOn()
72march_cubes.SetValue(0,260)
73march_cubes.Update()
74
75# Polydata mapper for the iso-surface
76isoMapper = vtk.vtkPolyDataMapper()
77isoMapper.SetInputConnection(march_cubes.GetOutputPort())
78isoMapper.ScalarVisibilityOff()
79
80# Actor for the iso surface
81isoActor = vtk.vtkActor()
82isoActor.SetMapper(isoMapper)
83isoActor.GetProperty().SetColor(1.0,1.0,1.0)
84
85# create iso metric renderer
86ren2.AddActor(isoActor)
87ren2.AddActor(txt2)

```

Figure 3 ISO metric rendering

```

90#####
91# View 3 for Volume+Iso mapping
92ren3 = vtk.vtkRenderer()
93txt3 = vtk.vtkTextActor()
94txt3.SetInput("ISO & Volume rendering!")
95txtprop=txt3.GetTextProperty()
96txtprop.SetFontFamilyToArial()
97txtprop.SetFontSize(18)
98txtprop.SetColor(1,1,1)
99
100# create a renderer for both volume as well as iso mapper
101ren3.AddVolume(volume1)
102ren3.AddActor(isoActor)
103ren3.AddActor(txt3)
104
105
106#####

```

Figure 4 Volume and ISO rendering

```

108# Render the scenes into 3 different viewports
109ren1.SetViewport(0, 0, 0.33, 1)    # viewport 1
110ren2.SetViewport(0.33, 0, 0.66, 1) # viewport 2
111ren3.SetViewport(0.66, 0, 1, 1)    # viewport 3
112
113# Rendering
114renWin = vtk.vtkRenderWindow()
115renWin.AddRenderer(ren1) # renders the volume mapper
116renWin.AddRenderer(ren2) # renders the iso mapper
117renWin.AddRenderer(ren3) # renders both volume and iso mapper
118
119# the render window size for 16:9 display
120renWin.SetSize(1920,1080)
121renWin.Render()
122
123# The window interactor
124iren = vtk.vtkRenderWindowInteractor()
125iren.SetRenderWindow(renWin)
126
127# syncing the viewports to first renderer
128camera = ren1.GetActiveCamera()
129ren1.ResetCamera()
130ren2.SetActiveCamera(camera)
131ren3.SetActiveCamera(camera)
132

```

Figure 5 Viewport division and rendering window

```

132
133# pixel intensities
134maximum_minimum = reader.GetOutput().GetScalarRange()
135print('the intensities are: ', maximum_minimum)
136# Dimension and voxel resolution
137dimension_voxel = reader.GetOutput(0)
138print(dimension_voxel)
139
140iren.Initialize()
141iren.Start()
142

```

Figure 6 Obtain the intensities, dimensions, size and voxel resolution

3. Output




Figure 7 Default view




Figure 8 Enlarged Right side view




Figure 9 Enlarged Side view




Figure 10 Enlarged Left side view

NOTE: The vtk version used for performing this program is **8.2.0**.

4. References

- 1) CT: <https://eclass.srv.ualberta.ca/mod/folder/view.php?id=3991376>
- 2) ISO: <https://eclass.srv.ualberta.ca/mod/folder/view.php?id=3991373>